

Appendix 3: Accuracy results


785B


Monovalent


DaKTP 005


DaKTP006

DTP134


DTP 135

CAIF SA  
 Dra. Bernarda Belay  
 Sub-Directora Técnica  
 M.O. 15.148

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

Appendix 4: Trendanalysis of 785B during the validation study

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


Appendix 5: Internal Standard: Comparison between the new method and the old method (on the basis of the most recent results)

	Type 1		Type 2		Type 3	
	New Method (DU/ml)	Old Method (DU/ml)	New Method (DU/ml)	Old Method (DU/ml)	New Method (DU/ml)	Old Method (DU/ml)
	37.2	37.1	6.9	7.9	25.6	27.2
	35.1	37.2	7.3	8.1	26.9	27.8
	35.8	36.4	7.8	8.2	32.3	26.5
	37.8	37.6	7.1	8.1	26.4	26.0
	35.1	37.9	6.5	7.7	25.4	25.8
	36.4	37.3	7.3	8.0	25.1	28.9
	38.1	36.7	7.7	8.3	26.1	25.2
	33.6	36.5	7.7	7.5	24.7	23.5
	33.4	36.8	7.9	8.0	27.0	23.9
	35.4	36.1	7.9	8.5	25.9	26.7
	35.8	34.2	7.7	7.9	24.0	25.4
	37.8	34.5	8.0	8.0	23.7	24.0
	37.2	37.0	8.5	7.7	25.9	24.5
	35.9	37.9	8.1	8.2	26.8	25.6
	32.0	36.2	7.2	8.0	19.5	26.1
	30.4	35.9	6.8	8.1	21.1	23.9
	40.4	36.3	7.7	8.3	26.2	25.6
	37.3	35.3	8.0	8.0	28.5	25.3
	39.2	36.8	8.9	8.0	31.1	25.6
	39.0	34.7	8.0	8.1	27.4	25.3
	44.0	36.6	8.5	8.3	29.6	25.8
Average	36.5	36.4	7.7	8.0	26.2	25.7
VC%	8.1	2.8	7.6	2.9	11.0	5.1
F-Test	0.00		0.00		0.00	
T-Test	0.91		0.06		0.49	

CAIF SA  
 Dra. Bernarda Belay  
 Co-Directora Técnica  
 M.A. 13-148  
 CAIF SA  
 Companhia Análises de  
 Investigaciones Farmacológicas S.A.  
 Dra. María Bernarda Belay  
 Analista  
 ONI 25378975

Report	RAP-20563 Version no.: 1
Retrospective validation of the sterility test	Page 1 of 15

**General information**

Old code N/A  
 Expiry period 25 years  
 Management Document management  
 Key words Sterility, retrospective, validation, growth promotion, antimicrobial activity

**Authorisation**

This document has been authorised by the document management system Quality On-line. The authorisation process was carried out according to procedure SOP-20091.

This document was compiled by ABC staff members L. Levels/J. Koopman

	Position	Name	Signature	Date
For consent	Document manager			

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

2 Identification ..... 4

2.1 Reference for method/SOP ..... 4

2.2 Recommendation for retrospective validation ..... 4

2.3 Review of methodology content and (current) validity ..... 4

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
RAP-20563

Version no.: 1

Retrospective validation of the sterility test

Page 3 of 15

## 1 Introduction

The sterility testing carried out within the LCB is described in a number of different documents that are closely related to one another and that form a single entity. For this reason all procedures have been incorporated in a single retrospective validation report.

### 1.1 Sterility testing

The sterility testing can be subdivided into direct inoculation (ANA-10157) and membrane filtration (ANA-10158). A number of controls should be complied with before the tests are declared valid:

- The growth promotion testing should demonstrate that the media are suitable for the growth of aerobic and anaerobic bacteria
- The negative controls should display no growth and be treated in exactly the same manner as the product being tested
- Conditions in the isolator should be controlled and should comply with regulations (ANA-10058)
- The product to be tested should be proven negative for antimicrobial activity.

Should these conditions not be met then the test is declared invalid and a repeat test must be carried out.

### 1.2 Growth promotion testing

The growth promotion test should show that the media used for sterility testing (Tryptone Soya Agar, Thioglycollate, Clausen) are suitable for the growth of aerobic and anaerobic bacteria and fungi. The following EP-prescribed test strains are used:

- Staphylococcus aureus
- Bacillus subtilis
- Pseudomonas aeruginosa
- Clostridium sporogenes
- Candida albicans
- Aspergillus niger

The test strains are prepared according to ANA-10068 and used to inoculate media at a concentration of 10-100 CFU. In practice, the media should comply with the following criteria:

- Visible bacterial growth within a maximum of 3 days
- Visible fungal growth within a maximum of 5 days
- Upper limit on count plates should be 150 CFU maximum
- For a lower limit less than 30 CFU per count plate, the test should be repeated to check for growth in the liquid medium

In practice, the sterility and growth promotion tests are carried out in parallel. Should the growth promotion test not meet all criteria, both tests should be repeated.

### 1.3 Antimicrobial activity

Before being tested for sterility, each product should be tested for antimicrobial activity in the applied media. In some cases, a product can have an effect on the growth promoting properties of the media. A growth-inhibiting effect renders the results of a sterility test unreliable.

Since each new product should be proven negative for antimicrobial activity before testing for sterility, this procedure can itself be considered part of product validation.

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RAP-20563

Version no.: 1

Retrospective validation of the sterility test

Page 4 of 15

Growth-inhibiting activity is determined at least *once* for each type of product (final lot/ final bulk) and repeated whenever the composition of the product or the composition (volume) of the relevant test media is changed.

## 2 Identification

### 2.1 Reference for method/SOP

The methods are described in the following protocols:

- EP §2.6.1. Sterility
- ANA 10058: Antimicrobial activity in presence of product (AAP)
- ANA 10157: Sterility testing
- ANA 10158: Sterility testing by membrane filtration
- ANA 10086: Growth promotion testing
- ANA 10065: Preparation of test suspensions

### 2.2 Recommendation for retrospective validation

The test method is limiting and compendial and has been carried out within the LCB since the early 1990's. It can be validated retrospectively. Since the validation is retrospective, there is no reference to a plan.

### 2.3 Review of methodology content and (current) validity

Compared to the procedure described in the EP, the applied method differs in three ways:


- For practical reasons, a distinction is made between  $<2$  ml and  $\geq 2$  ml samples. This is done so that the well-established sampling regimen within the LCB need not be changed, enabling the use of old validation data for new results. In addition, a larger sample volume aids visibility of the samples during the procedures in the isolator. Use of a larger sample volume has no effect on the test result.
- Clausen medium is used instead of TSB when Merthiolate is present in the product.
- A maximum of 150 CFU is maintained when using the test strains. Given the range in results obtained, this limit is more easily maintained than the EP recommendation of 10-100 CFU. A 20-50% range in results is acceptable for microbiological culture methods.

The EP does allow some deviation from the recommended procedure, certainly if this results in direct improvements and therefore in (further) improvement in the reliability of the results. N.B. Authorisation of the relevant ANA's implies inherent approval of this deviation from EP recommendations.

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RAP-20563

Version no.: 1

Retrospective validation of the sterility test

Page 5 of 15

## 2.4 Requirements and criteria

Since the sterility testing within the LCB is described in a number of different documents, all these procedures are to be included in a single validation. The test analysis follows the procedures described in the European Pharmacopoeia, meaning that the tests are compendial, that they are limiting and that they can be validated retrospectively. The aim of a retrospective test is to guarantee the continuity of a testing method and thereby its consistency. In addition, data are presented that show that the analysis also complies with the ICH parameters for a limiting test regarding specificity and detection limit

## 2.5 Abbreviations

LCB	- Laboratory for the Control of Biological products
ICH	- International Conference on Harmonisation
EP	- European Pharmacopoeia
TSB	- Tryptone Soya Broth
Thio	- Thioglycollate
ATCC	- American Type Culture Collection
AAP	- Antimicrobial activity in presence of product
CFU	- colony-forming unit
m.o.	- microorganisms
RIA	- request for amendment

## 3 Results

### 3.1 Specificity

The specificity of the analysis can be separated into three components:


- specificity of the applied media
- specificity/identity of the test strains
- the presence during analysis of antimicrobial activity from the vaccine itself, which might influence the specificity of the media/m.o.

N.B. The specificity (composition) of the applied media and the test strains is registered in the European Pharmacopoeia.


#### 3.1.1 Media

The applied media have a standard composition in accordance with the EP. The specificity of the media used in the sterility testing implies that the media are chosen to allow the isolation of a broad range of aerobic and anaerobic m.o. The descriptions of the applied media demonstrate that this is the case:

- Tryptone Soya Broth (TSB) – a general medium containing two peptones that promotes the growth of a wide range of microorganisms (m.o.). Primarily suitable for the culture of aerobic m.o., but also suitable for the culture of anaerobic m.o. when used under anaerobic conditions.
- Thioglycollate (Thio) – primarily suitable for the culture of anaerobic m.o., but aerobic m.o. will also grow in this medium.
- Clausen – a medium containing neutralising components and essential metals for the recovery of anaerobic spore-forming bacteria; specially developed for products containing Merthiolate.

Report

RAP-20563

Version no.: 1

Retrospective validation of the sterility test

Page 6 of 15

## 3.1.2 Antimicrobial activity

Before a product can be tested for sterility, it should first be shown to contain no antimicrobial activity. The products shown in table 1 are some of those which have undergone testing for antimicrobial activity and for which the test specificity has been analysed.

Table 1: results of AAP testing

Culture temp.	20 - 25 °C										30 - 35°C			
	A					B					C	D	E	F
Strain	TSB		Claus		TSB		Thio		Claus		Thio			
Culture days	3	5	3	5	3	5	3	5	3	5	3			
DTP	0		*	*	0		*	*	*	*	C	0	0	0
DKTP	0		*	*	0		*	*	*	*	0	0	0	0
BMR	0		*	*	0		*	*	*	*	0	0	0	0
Polio	0		*	*	0		*	*	*	*	0	0	0	0
BCG	0		*	*	0		*	*	*	*	0	0	0	0
Acell. pertussis		0	*	*	0		*	*	*	*	0	0	0	0
Tetanus	S		0		S		S	S	S	S	0	0	0	C

A - Candida albicans  
 B - Aspergillus niger  
 C - Staphylococcus aureus  
 D - Clostridium sporogenes  
 E - Pseudomonas aeruginosa  
 F - Bacillus subtilis  
 O - no inhibition  
 M - growth inhibited  
 S - growth strongly inhibited  
 \* - N/A

## 3.1.3 Test strains (control organisms)

Six different test strains are used as positive controls within the AAP and growth promotion tests. Use of these strains guarantees test-specificity for a wide range of m.o.

All six strains are EP-approved and have an ATCC number (see table 2).

Before routine use in analysis, test suspensions are made according to ANA-10065. As an extra control, these suspensions are checked for their identity to verify that the correct strains are in use. See appendices 1 to 6 for the results obtained.

Table 2: Control organisms for sterility testing

Microorganism	ATCC no.	Description
Candida albicans	10231	Yeast
Aspergillus niger	16404	Fungus
Staphylococcus aureus	6538	Aerobic gram-positive bacteria
Clostridium sporogenes	19404	Anaerobic spore-forming, gram-positive bacteria
Pseudomonas aeruginosa	9027	Aerobic gram-negative bacteria
Bacillus subtilis	6633	Aerobic spore-forming, gram-positive bacteria


## 3.2 Detection limit

For a limiting test, the detection limit should be validated. During sterility testing, the detection limit of each separate test is determined using a growth promotion test. This demonstrates that the applied media are able to detect bacteria in the range of 10 - 150 CFU. Should the growth promotion test not fall within this range, then the whole test is declared invalid and a repeat test carried out using a different batch of medium.

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3.3 Continuity

Continuity can be demonstrated using results obtained with the test strains in growth promotion tests (in the period Jan. 2002 - Sept. 2003). At the start of the growth promotion test, two count plates are inoculated in order to confirm the number of microbes in the test suspensions. The average count from the count plates is shown in figures 1 to 6 for each strain.

Figure 1: Counts for Staphylococcus aureus


Figure 2: Counts for Bacillus subtilis


Figure 3: Counts for Candida albicans


Figure 4: Counts for Clostridium sporogenes


Figure 5: Counts for Pseudomonas aeruginosa



Figure 6: Counts for Aspergillus niger

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RAP-20563

Version no.: 1

Retrospective validation of the sterility test

Page 8 of 15

## 4 Discussion

### 4.1 Specificity

All six of the test strains demonstrated abundant growth in the presence of products with preservatives (DKTP, DTP, Polio, Acellular pertussis) and products without preservatives (BCG, BMR), in all tests. This growth was visually no different to the incorporated controls both in TSB at 20-25°C and in Thio at 30-35°C. The controls consist of the test strains grown in TSB and Thio in the absence of product.

In the Tetanus vaccine, however, the presence of Merthiolate leads to growth inhibition in two of the six test strains: *Candida albicans* and *Aspergillus niger*. Clausen medium can be used to replace TSB at 20-25°C to solve this problem since it has a neutralising effect on Merthiolate. Using this medium for sterility tests of Tetanus vaccines solves the problem of growth inhibition in *Candida albicans*; the growth of *Aspergillus niger* is however still inhibited. Despite further analysis of other media and the addition of Na-Thioglycollate to TSB medium at 20-25°C in an attempt to neutralize Merthiolate, a solution has not been found. For this reason, a request for amendment (RfA) has been instigated that states that the growth of *Aspergillus niger* may be left out of the results of AAP testing in the presence of Merthiolate.

Since this relates to only one of the six control organisms, and this microorganism is known to have tricky growth properties, we have no doubts as to the validity of the test. The addition of Merthiolate as a preservative for the Tetanus vaccine certainly reduces the chances of growth and presence of contaminating m.o. in the product. In an earlier version of the AAP test where Clausen medium was also applied, this problem did not present itself since at that time *Aspergillus niger* was not a test strain required by the EP.

### 4.2 Detection limit and continuity

For each separate test, the detection limit is determined in the growth promotion tests of the applied media. The use of the six test strains forms a basis for evidence of the continuity of the analysis. Figures 1 to 6 demonstrate the stability of the test strains during the period shown. In the graphs of the *Bacillus subtilis* trend line (one point) and the *Pseudomonas aeruginosa* trend line (two points), there are three points that deviate from the rest. In these tests, a double amount of test suspension was added by accident; they were however declared valid since the remaining controls were satisfactory.

The trend line for *Clostridium sporogenes*, cultured in an anaerobic medium, is more scattered. This strain spreads unevenly on the plate, making it difficult to count/estimate and therefore giving a wide range of results.

## 5 Conclusion

Based on the results of this retrospective validation of testing for the presence of contaminating microorganisms, as carried out within the LCB since the early 1990's, we can conclude that we have shown continuity of the analysis and have also shown a guaranteed consistency.

## 6 Appendices

Appendices 1 - 6: Identity of test strains

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**Identification report carried out by**

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Your reference number: order number 304334

Our reference number: det. 301, 1999

**Description of the identification method:**

-media used: Czapek agar, Malt Extract agar (Oxoid)

-temperature of incubation: 25 °C

-incubation: 7-14 days in darkness

-scientist who carried out the identification : E.S. Hoekstra


- (main) reference/description used for identification:

Raper, K.B. and Fennell, D.I.(1965) The Genus *Aspergillus*. Baltimore, Maryland; Williams and Wilkins**Identification criteria:**


Identification was based on colony morphology and microscopical examination using light microscopy.

Postal address: P.O. Box 273, 3740 AG BAARN, The Netherlands  
 Visiting address: Oosterstraat 1, 3742 SK BAARN, The Netherlands

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Page 15 of 15

**Short description of the isolate.**

Colonies on Czapek agar attaining a diameter of 5-6 cm diameter in 10 days at 25 °C, consisting of a blackish-brown felt due to sporulating heads.

Colonies on Malt extract agar, attaining a diameter of 4,5-5 cm after 10 days, consisting of a blackish-brown felt due to sporulating heads and also whitish aerial mycelium could be observed.

**Microscopical observation:** Conidial heads, biserial, blackish-brown, radiate and splitting into columns with age. Conidiophore stipe, smooth-walled, hyaline often brownish towards the vesicle. Vesicles (sub)globose, mostly 30-50 µm in diam. Conidia (sub)globose, blackish, commonly 5 µm in diam, spiny,

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Investigaciones Farmacéuticas S.A.