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

RAP-21890

Versie: 1

Determination of host cell proteins during production of monovalent pools IPV

Blad 2 van 17

Table of contents

1.	Summary	3
2.	Introduction	4
3.	Materials and methods	5
3.1.	Production of HCP mix and polyclonal anti HCP-mix antibodies	5
3.2.	Sample materials	5
3.3.	Dot Blotting	6
3.3.1.	Reagents	6
3.3.2.	Sample dilutions	6
3.3.3.	Preparation of calibration standards	7
3.3.4.	Dot-blot procedure	7
3.3.5.	Data analysis	7
3.4.	SDS page and Western Blotting (used as control test)	8
3.4.1.	Reagents	8
3.4.2.	SDS page and Western blot procedure	8
4.	Results and Discussion	9
4.1.	Qualification of dot blot test	9
4.1.1.	Specificity of polyclonal antibodies	9
4.1.2.	Linearity and range	10
4.1.3.	Reproducibility and accuracy	14
4.1.4.	Conclusion qualification of dot blot analysis	14
4.2.	HCP determination in IPV production fraction 2.1, 3.1 and 5.1.	14
4.2.1.	IPV production fraction 2.1	14
4.2.2.	IPV production fraction 3.1	14
4.2.3.	IPV production fraction 5.1	15
4.3.	Clearance of HCP during production	15
5.	Conclusions	17
6.	Archive location raw data	17

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RAP-21890

Versie: 1

Determination of host cell proteins during production of
monovalent pools IPV

Blad 3 van 17

1. Summary

A dot blot method has been developed for quantitative determination of host cell proteins (HCP,s) during the production of inactivated poliomyelitis vaccine (IPV).

The method has been qualified as follows:

- quantification range of 0.19-3.8 µg/mL
- accuracy is 80-120%
- Reproducibility as CV <30%

HCPs were measured in poliomyelitis virus culture samples of 9 production fractions 2.1 (after clarification and concentration), 9 production fractions 3.1 (after size exclusion chromatography) and 9 production fractions 5.1 (after ion exchange chromatography and sterile filtration). Three batches of each type of poliomyelitis vaccine were used.

On average production fraction 2.1 samples contain 13 mg/mL HCP, production fraction 3.1 samples contain 0.27 mg/mL HCP, while samples of production fraction 5.1 contain less than 0.5 µg/mL HCP. Since the volume in fraction 5.1 is comparable to the volume in fraction 6.1, we conclude that monovalent pools of IPV contain less than 0.5 µg/mL HCP. No significant differences could be observed for the three different types of poliovirus productions.


Clearance of HCP after size exclusion chromatography and ion exchange chromatography was higher than 99.9% in all tested batches.

It can be concluded that during IPV production HCP is cleared efficiently by the combination of size exclusion chromatography and ion exchange chromatography. The removal of HCP is consistent for all three types of poliovirus.

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RAP-21890

Versie: 1

Determination of host cell proteins during production of
monovalent pools IPV

Blad 4 van 17

2. Introduction

During poliovirus cultivation the vero host cells lyse. Proteins are released into the medium and cell fragments detach from the micro carriers. Cell fragments and host cell proteins are removed during down stream processing (DSP). DSP consists of filtration and chromatographic purification steps (Table 1). Traces of host cell proteins (HCP's) however will remain in the monovalent pools. We investigate the amount of residual host cell proteins and the effect of the filtration and purification steps on the removal of the host cell proteins.

In order to do so it was decided to use the technique of dot blot. Shortly, dot blotting consists of the following steps:

1. Sample or standard materials are spotted on a nitrocellulose membrane
2. The membrane is incubated with anti-HCP antibody
3. The membrane is washed to remove non-specifically bound materials
4. The membrane is incubated with second antibody labeled with horseradish peroxidase (HRP)
5. The membrane is washed to remove non-specifically bound materials
6. As an HRP signal reagent a chemiluminogenic reaction based on GZ-11 is used (Quantitative western blotting using a GZ-11 based chemiluminogenic signal reagent. G Zomer *et al* in Proceedings of the 14th International Symposium on Bioluminescence and Chemiluminescence. Ed AA Szalay, PJ Hill LJ Kricka PE Stanley 2006: 171-4)
7. After addition of the signal reagent the glowing membrane is imaged and the resulting image of the blots is quantified

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RAP-21890

Versie: 1

Determination of host cell proteins during production of
monovalent pools IPV

Blad 5 van 17

Table 1. Production stages of IPV monovalent pools

Code	Stage
Fraction 0.1	After cultivation
Fraction 1.1	After clarification (0.22 µm)
Fraction 2.1	After concentration 100 kD
Fraction 3.1	After size exclusion chromatography
Fraction 4.1	After ion exchange chromatography
Fraction 5.1	After sterile filtration (0.22 µm)
Fraction 6.1	After inactivation

3. Materials and methods

3.1. Production of HCP mix and polyclonal anti HCP-mix antibodies

The polyclonal rabbit anti-HCP serum was obtained by a cascade immunization of rabbits with a HCP mix from Vero cells. Three days before immunization blood was obtained from the rabbits. This pooled serum serves as pre-immune rabbit serum (NKS)

- The HCP mix, a blank culture of Vero cells is 3 times frozen and thawed to induce cell lysis. The culture supernatant was pooled and filtered. The pore size of the filters was comparable with these used during the IPV Vero production process. The HCP mix was aliquot, rapid frozen using dry ice and stored at -70°C until use.
- Rabbit anti-HCP serum: 3 rabbits were immunized with the HCP-mix and received a booster with the same HCP mix after 28 days. On day 42 blood samples were collected. Using the antibodies as an affinity column the immunodominant proteins have been removed from the HCP-mix. The treated HCP-mix is used to immunize the 3 rabbits again on day 56. The same procedure was performed with serum of day 77 to immunize at day 85. The rabbits were given a booster on day 121 followed by collection blood on day 135. The pooled serum of day 135 is used in this study as anti-HCP serum.

3.2. Sample materials

The poliomyelitis virus culture samples that were studied are presented in table 1. Fraction 2.1 is the fraction obtained after clarification (filtration through a three-step in-line filter unit containing 75 µm, 0.45 µm and 0.22 µm filters) and concentration (over a 100kD cut off filter), fraction 3.1 is the fraction obtained after purification (using gel filtration and anion exchange chromatography), fraction 5.1 is the fraction after addition of medium 199 and sterile filtration). Fraction 5.1 was used instead of 6.1 because the latter is formaldehyde treated. The reactivity of the antiserum may be different with formaldehyde treated material. Because the HCP content will not change from 5.1 to 6.1, it was decided to test upstream, i.e. in fraction 5.1.

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
RAP-21890

Versie: 1

**Determination of host cell proteins during production of
monovalent pools IPV**

Blad 6 van 17

Table 2: IPV sample materials used during this study

batch number	Monovalent pool IPV	Polio type	Production fraction	Production fraction	Production fraction
PV07	131	1	2.1	3.1	5.1
PV07	132	1	2.1	3.1	5.1
PV07	133	1	2.1	3.1	5.1
PV07	216	2	2.1	3.1	5.1
PV07	217	2	2.1	3.1	5.1
PV07	218	2	2.1	3.1	5.1
PV07	325	3	2.1	3.1	5.1
PV07	326	3	2.1	3.1	5.1
PV07	327	3	2.1	3.1	5.1

3.3. Dot Blotting

3.3.1. Reagents

- Wash buffer: Phosphate buffered saline 10mM, pH 7.3 (PBS) containing 0.1% Tween20, pH 7.3 (PBST)
- Assay buffer: Wash buffer containing 0.5% (w/v) of Protifar
- Polyclonal anti-HCP (see 2.1)
- 2ndAb.HRP (monoclonal anti rabbit IgG (γ -chain specific) HRP conjugate

3.3.2. Sample dilutions

- IPV production fraction 2.1 samples (see table 1) are initially diluted about 5000x using weight of samples and diluents. From this stock, four independent dilutions are prepared using calibrated pipettes as follows (table 3):


Table 3: preparation of 2.1 sample dilutions

μ l of sample	μ l of PBS
50	450
150	350
250	250
350	150

- IPV production fraction 3.1 samples (see table 2) are initially diluted about 100x using weight of samples and diluents. From this stock, four independent dilutions are prepared using calibrated pipettes as follows (table 3).
- IPV production fraction 5.1 samples are used undiluted, and after pretreatment with 0.015% hydrogen peroxide (final volume) in order to pre-oxidize possibly interfering anti-oxidants that were present in the samples

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
RAP-21890

Versie: 1

**Determination of host cell proteins during production of
monovalent pools IPV**

Blad 7 van 17

3.3.3. Preparation of calibration standards

- HCP-mix (concentration 0.76 mg/mL, determined by a biocinchonic acid assay) is initially diluted about 100x using weights of standard and diluent. From this stock, seven independent dilutions and a blank are prepared using calibrated pipettes as follows (table 4).

Table 4: preparation of HCP calibration standards

µl of stock	µl of PBS
500	500
400	600
300	700
200	800
100	900
50	950
25	975
0	1000

3.3.4. Dot-blot procedure

- Suitably diluted (in PBS) samples (see 2.3.1 and 2.3.2) or standards (see 2.3.3, 50µL) were applied to a nitrocellulose membrane using an ELIFA Unit.
- After washing with 200 µL of PBS buffer, the membrane was incubated with polyclonal anti-HCP (see 2.1, 1:300 in assay buffer) during one hour with gentle shaking.
- After washing with wash buffer, the membrane was incubated with 2ndAb.HRP (1:100,000 in assay buffer) during one hour with gentle shaking.
- The membrane was washed (4 x) with wash buffer.
- Signal reagent was prepared by adding 4-phenylphenol (5 mg/mL, 40µL) and GZ-11 (7 mg/mL, 50 µL) to 5 mL of PBS containing 0.1 mM EDTA and 0.1 % Triton-X-100. To the signal reagent was added 5 mL of Na-acetate buffer, 160 mM, pH 6, containing 0.015% of hydrogen peroxide. The resulting mix was added to the membrane, and after incubating for 5 minutes, decanted. The glowing membrane was measured in a chemiluminescent imager. Resulting images were stored as TIF-files.

3.3.5. Data analysis

- Image files are loaded into Image J software.
- The image is inverted giving dark dots against a light background
- Rectangular row selections are made and converted to profiles
- From the profiles peak areas are measured
- Peak areas are exported to MS-Excel
- MS-Excel routines are used for regression, 4 parameter fits, and further calculations. Sample concentrations were calculated from 5 out of the 8 available data points per sample. The three data points that differed most from the average were excluded from the calculation. Because of the low HCP-concentrations in 5.1 samples, a limited number of calibration standards ranging from 0.16-0.64 µg/mL were used. In this concentration range a linear calibration curve could be used.

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RAP-21890

Versie: 1

Determination of host cell proteins during production of
monovalent pools IPV

Blad 8 van 17

3.4. SDS page and Western Blotting (used as control test)

3.4.1. Reagents

- Electrophoresis buffer: Tris 25mM, Glycine, 192 mM, SDS, 0.1% ((w/v), pH 8.3.
- SDS-gels: 12% Tris-HCl precast
- Molecular weight markers ranging from 31-200kD
- Western Blot buffer Tris 25mM, Glycine 192 mM, Methanol (20%, v/v), pH 8.3
- Nitrocellulose membrane
- Wash buffer: Phosphate buffered saline 10mM, pH 7.3 (PBS) containing 0.1% Tween20, pH 7.3 (PBST)
- Assay buffer: Wash buffer containing 0.5% (w/v) of Protifar
- Polyclonal anti-HCP (see 2.1)
- 2ndAb.HRP (monoclonal anti rabbit IgG(γ-chain specific) HRP conjugate


3.4.2. SDS page and Western blot procedure

- A suitable amount of sample was first precipitated using a 72% w/w trichloroacetate (TCA)-solution. The final concentration of TCA after mixing with sample is 7.2%.
- The samples were thoroughly vortexed followed by centrifugation
- The supernatant was carefully removed and 15 µL of reducing buffer was added.
- The samples were thoroughly vortexed and placed in a thermostat at 100°C during 5 minutes. After this the samples were centrifuged for a few seconds, only to let the liquid in the top of the sample container to be spun down.
- The slots of the precast gel were loaded with 15 µL of sample. Molecular weight markers used as a reference were loaded with 7 µL .
- After electrophoresis the gel was washed with water and either coloured with a coomassie blue staining or the western blot procedure as follows was started.
- The gel was washed again with blot buffer
- Via a cold tank transfer in blotting buffer a Western Blotting was performed
- After washing with µLPBS buffer, the membrane was incubated with polyclonal anti-HCP (see 2.1, 1:100 in assay buffer) during one hour with gentle shaking.
- After washing with wash buffer, the membrane was incubated with 2ndAb.HRP (1:100,000 in assay buffer) during one hour with gentle shaking.
- The membrane was washed (4 x) with wash buffer.
- Signal reagent was prepared by adding 4-phenylphenol (5 mg/mL, 40µL) and GZ-11 (7 mg/mL, 50 µL) to 5 mL of PBS containing 0.1 mM EDTA and 0.1 % Triton-X-100. To the signal reagent was added 5 mL of Na-acetate buffer, 160 mM, pH 6, containing 0.015% of hydrogen peroxide. The resulting mix was added to the membrane, and after incubating for 5 minutes, decanted. The glowing membrane was measured in a chemiluminescent imager.
- The resulting image was saved as a TIF-file and analyzed using Image-J software.

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RAP-21890

Versie: 1

Determination of host cell proteins during production of
monovalent pools IPV

Blad 9 van 17

4. Results and Discussion

4.1. Qualification of dot blot test

4.1.1. Specificity of polyclonal antibodies

From Western blots (figure 1) and SDS-PAGE gel (figure 2) it is clear that the polyclonal anti-HCP antibodies recognize vero host-cell proteins. In the western blot the HCP mix shows two main bands at about 35-38 kD together with minor bands e.g. at about 50 kD after staining with polyclonal anti-HCP antibodies. The production fraction 2.1 (figure 1, lane 2) and production fraction 3.1 (figure 1, lane 3) samples show almost super imposable bands, most of which are also present in the HCP-mix (figure 1, lane 4). Differences can probably be attributed to matrix effects (the HCP-mix was isolated from vero cells which had not been exposed to virus). The production fraction 5.1 sample (figure 1, lane 5) shows no HCP-related bands. The inverted bands indicated by VP1-3 are virus related proteins (see also figure 2), two of which are also visible in the 3,1 fraction (figure 1, lane 3). Lane 6 in figure 1 (normal rabbit serum) shows no response, indicating the absence of HCP in normal rabbit serum.



Figure 1: Western blot stained with HCP polyclonal antibodies.

Lane 1 contains mw markers, lane 2 is a sample of production fraction 2.1 (PV07.218) diluted 250 times, lane 3 contains a sample of production fraction 3.1 (PV07.218) diluted 10 times, lane 4 contains the HCP-mix diluted 10 times, lane 5 contains an undiluted sample of production fraction 5.1 (PV07.218). Lane 6 is empty and lane 7 contains normal rabbit serum as a negative control. Lane 8 is empty

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
RAP-21890

Versie: 1

Determination of host cell proteins during production of
monovalent pools IPV

Blad 10 van 17


Figure 2: SDS-Page electrophoresis gel of Figure 1 after residual protein staining. Lane 1 protein markers, lane 2: empty, lane 3: magic marker (corresponds with lane 1 of figure 1), lanes 4-8 correspond with lanes 2-6 of figure 1. lane 4 is a sample of production fraction 2.1 (PV07.218) diluted 250 times, lane 5 contains a sample of production fraction 3.1 (PV07.218) diluted 10 times, lane 6 contains the HCP-mix diluted 10 times, lane 7 contains an undiluted sample of production fraction 5.1 (PV07.218). Lane 8 is empty while lane 9 contains normal rabbit serum as a negative control. Lane 10 is empty

When comparing figures 1 and 2 it is clear that during the purification process HCP's are cleared from the product: while virus proteins VP1-3 are faintly visible in production fraction sample 2.1 (Figure 2, lane 4), these proteins become increasingly abundant in production fraction sample 3.1 (figure 2, lane 5) and production fraction sample 5.1 (figure 2, lane 7) at the expense of host cell proteins.

4.1.2. Linearity and range

The dot blot in figure 3 shows dilution ranges of HCP mix and 2.1 and 3.1 IPV production fraction samples stained with HCP polyclonal antibodies.
The dot blot shows visible blots in a HCP mix range of 3.8-0.19 µg/mL.

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RAP-21890

Versie: 1

Determination of host cell proteins during production of monovalent pools IPV

Blad 11 van 17


Figure 3: Dot blot stained with HCP polyclonal antibodies.

The top two rows are calibration standards of HCP mix in duplicate (7 dilutions from 3.8-0.19 µg/mL plus blank, from left to right). The third and fourth rows are 8 dilutions in duplicate (48,000-6,000, left to right) of a 2.1 IPV production fraction sample (PV-07.218), while the fifth and sixth rows are 8 dilutions in duplicate (2,000-250, left to right) of a 3.1 IPV production fraction sample (PV07.218).

Using ImageJ software the intensity of the dots can be quantified. In this example the dots in each row are converted to profile data as indicated in figure 4:




Figure 4: profile data of dot blot from figure 3.

From these profiles the areas of the peaks can be measured. The top two rows are calibration standards in duplicate (7 dilutions from 3.8-0.19 µg/mL plus blank, from left to right). The third and fourth rows are 8 dilutions in duplicate (48,000-6,000, left to right) of a 2.1 production fraction sample (PV-07.218), while the fifth and sixth rows are 8 dilutions in duplicate (2,000-250, left to right) of a production fraction 3.1 sample (PV07.218).

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RAP-21890

Versie: 1

Determination of host cell proteins during production of
monovalent pools IPV

Blad 12 van 17

From the profile data peak areas can be measured for each peak. The data are shown in table 5:

Table 5: peak areas (in arbitrary units (AU)) as obtained from the profile data of figure 4. The top two rows are calibration standards in duplicate (7 dilutions from 3.8-0.19 µg/mL plus blank, from left to right). The third and fourth rows are 8 dilutions in duplicate (48,000-6,000, left to right) of a 2.1 production fraction sample (PV-07.218), while the fifth and sixth rows are 8 dilutions in duplicate (2,000-250, left to right) of a production fraction 3.1 sample (PV07.218).


Cal.	4296	3255	2135	1464	929	483	231	0
Cal.	3852	2792	1969	1237	664	414	242	0
2.1	342	546	710	898	1042	1281	1418	1529
2.1	251	487	834	920	1011	1251	1497	1544
3.1	193	381	405	518	677	777	815	922
3.1	215	303	353	544	588	648	793	811


A calibration curve can be constructed by plotting the peak areas of the calibration standards against their concentrations (see table 6 for HCP calibration concentrations and sample dilution factors).

Table 6: Concentrations of HCP calibration standards in duplicate (rows 1-2), together with dilution factors of samples PV07.218 production fraction 2.1 in duplicate (rows 3-4) and production fraction PV07.218 3.1 in duplicate (rows 5-6), respectively

3.8	3.04	2.28	1.52	0.76	0.38	0.19	0	HCP mix dilutions µg/mL
3.8	3.04	2.28	1.52	0.76	0.38	0.19	0	HCP mix dilutions µg/mL
48228	24114	16076	12057	9646	8038	6890	6028	Fraction 2.1 dilutions
48228	24114	16076	12057	9646	8038	6890	6028	Fraction 2.1 dilutions
2034	1017	678	508	407	339	291	254	Fraction 3.1 dilutions
2034	1017	678	508	407	339	291	254	Fraction 3.1 dilutions

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Determination of host cell proteins during production of monovalent pools IPV

The calibration curve is shown in figure 5:


Figure 5: HCP calibration curves using linear regression and 4 parameter fit


The linear range of the calibration line extends from 0.19 to 3.8 µg/mL. From this calibration curves sample concentrations can be calculated from their responses (see table 5), taking into account their dilution factors (see table 6). In this example the concentration HCP in each of the two samples can be calculated from 16 data points (eight dilutions in duplicate) using either linear or four parameter calibration. Results are shown in table 7:

Table 7: sample concentration using linear or four parameter fit using all data (n=16) or data without outliers (n=10)

n=16	linear calibration		4 parameter calibration	
	218 fraction 2.1	218 fraction 3.1	218 fraction 2.1	218 fraction 3.1
sample				
conc (µg/mL)	11741	310	11546	259
%CV	21%	31%	15%	18%

n=10	linear calibration		4 parameter calibration	
	218 fraction 2.1	218 fraction 3.1	218 fraction 2.1	218 fraction 3.1
sample				
conc (µg/mL)	10366	262	10834	252
%CV	8%	9%	7%	6%

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RAP-21890

Versie: 1

Determination of host cell proteins during production of
monovalent pools IPV

Blad 14 van 17

4.1.3. Reproducibility and accuracy

Reproducibility was deduced from repeating the analysis (data not shown). Variation coefficient is <30%. Accuracy was assessed from back-calculated concentrations of the standard line (data not shown). The accuracy is 80-120%.

4.1.4. Conclusion qualification of dot blot analysis

The dot blot analysis provides a test method with adequate sensitivity (0.19 µg/mL), linear range (0.19-3.8 µg/mL), reproducibility (CV<30%), and accuracy (80-120%) to allow quantification of HCP's in production samples 2.1, 3.1, and 5.1.

4.2. HCP determination in IPV production fraction 2.1, 3.1 and 5.1.

4.2.1. IPV production fraction 2.1

HCP concentrations in 9 IPV production fraction 2.1 samples are determined and presented in Table 8. The HCP concentrations are in the order of 10-20 mg/mL. No significant differences in HCP amount are observed between the different types of poliovirus. Protein concentrations as determined using a biocinchonic acid assay are in the same range (data not shown) Fraction 2.1 contains thus almost exclusively HCP.

Table 8: Results of HCP dot blot assay in IPV production fraction 2.1 samples.
The average of two dot blot assays is calculated

batch number	Polio type	HCP conc. in µg/mL (average)
PV07-131	1	7475*
PV07-132	1	10039
PV07-133	1	16449
PV07-216	2	13829
PV07-217	2	12185
PV07-218	2	11475
PV07-325	3	11504
PV07-326	3	13980
PV07-327	3	18677

* Sample is only tested once. An average could not be determined

4.2.2. IPV production fraction 3.1

HCP concentrations in 9 IPV production fraction 3.1 samples are determined and presented in Table 9. HCP concentrations range from about 0.19 – 0.35 mg/mL. No significant differences in HCP amount are observed between the different types of poliovirus. From the reduction of HCP concentration in fraction 3.1 compared to the HCP concentration in fraction 2.1, the conclusion can be drawn that size exclusion chromatography removes HCP (see also table 10). However still a remaining amount of HCP is present in fraction 3.1.

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RAP-21890

Versie: 1

Determination of host cell proteins during production of
monovalent pools IPV

Blad 15 van 17

Table 9: Results of the HCP dot blot assay in IPV production fraction 3.1 samples.
The average of two dot blot assays is calculated

batch number	Polio type	HCP conc. in µg/mL (average)
PV07-131	1	301
PV07-132	1	198
PV07-133	1	283
PV07-216	2	289
PV07-217	2	276
PV07-218	2	187
PV07-325	3	258
PV07-326	3	244
PV07-327	3	357

4.2.3. IPV production fraction 5.1

HCP concentrations in 9 IPV production fraction 5.1 samples are determined and presented in Table 10. HCP-concentrations in 5.1 samples are below 0.5 µg/mL in all 9 tested samples. HCP concentration is independent of poliovirus type. The results show that the ion exchange chromatography contributes to a good clearance of HCP (see also table 11).

Table 10: Results of the HCP dot blot assay in IPV production fraction 5.1 samples.

batch number	Polio type	HCP conc. µg/mL
PV07-131	1	0.11
PV07-132	1	0.07
PV07-133	1	0.10
PV07-216	2	0.32
PV07-217	2	0.14
PV07-218	2	0.44
PV07-325	3	0.14
PV07-326	3	0.12
PV07-327	3	0.19

4.3. Clearance of HCP during production


Table 11 gives an overview of the measured HCP-concentrations in IPV production fractions 2.1, 3.1, and 5.1. From this table it can be concluded that during production HCP is cleared efficiently by the combination of size exclusion chromatography and ion exchange chromatography. The removal of HCP is consistent for all three types of poliovirus.

CAIF SA
Dra. Bernarda Belay
Co-Directora Técnica
M.N. 15.148
Compañía Argentina de
Investigaciones Farmacéuticas S.A.
Dra. ~~Maria Bernarda Belay~~
Apoderada
DNI 29378925

NVI
nederlands vaccin instituut

Rapport

RAP-21890

Versie: 1

Determination of host cell proteins during production of
monovalent pools IPV

Blad 16 van 17

Table 11: Overview of HCP-concentrations in production fraction 2.1, 3.1, and 5.1

Batch number	Polio type	HCP concentrations (µg/mL)		
		Fraction 2.1	Fraction 3.1	Fraction 5.1
PV07-131	1	7475	301	0.11
PV07-132	1	10039	198	0.07
PV07-133	1	16449	283	0.10
PV07-216	2	13829	289	0.32
PV07-217	2	12185	276	0.14
PV07-218	2	11475	187	0.44
PV07-325	3	11504	258	0.14
PV07-326	3	13980	244	0.12
PV07-327	3	18677	357	0.19
average		12800	266	0.18
SD		2410	52	0.12

D antigen concentration are determined as in process controles during production. HCP clearance can be determined when the HCP amount is calculated in relation to the D antigen amount. Results for the process step going from 3.1 to 5.1 are shown in table 12.

Table 12: HCP clearance after ion exchange chromatography.
Clearance is calculated as $(1 - (\text{HCP}_{5.1}/\text{D-antigen}_{5.1} - \text{HCP}_{3.1}/\text{D-antigen}_{3.1})) * 100\%$

	3.1 D antigen /mL	3.1 HCP (µg /mL)	5.1 D antigen /mL	5.1 HCP (µg /mL)	clearance 3.1->5.1
PV07-131	7836	301	2548	0.11	99.9%
PV07-132	7690	198	2636	0.07	99.9%
PV07-133	na	283	na	0.10	na
PV07-216	2838	289	849	0.32	99.6%
PV07-217	2536	276	1090	0.14	99.9%
PV07-218	1745	187	937	0.44	99.6%
PV07-325	5647	258	1440	0.14	99.8%
PV07-326	5127	244	1243	0.12	99.8%
PV07-327	4416	357	1302	0.19	99.8%

na=not available

HCP clearance after ion exchange chromatography was higher than 99.6% in all tested batches. No significant differences in clearance could be observed for the three different types of poliovirus productions.

HCP clearance after two purification steps (after size exclusion chromatography and ion exchange chromatography) can be calculated in the same way using the data of table 13. As can be seen from table 13 HCP is removed efficiently >99.9% after size exclusion and ion exchange chromatography. No significant differences in clearance could be observed for the three different types of poliovirus productions.

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