


Customer Name \_\_\_\_\_ Customer No. \_\_\_\_\_ Purchase Order No. \_\_\_\_\_ Fax \_\_\_\_\_

Billing Address \_\_\_\_\_ Patient Name \_\_\_\_\_

Shipping Address \_\_\_\_\_ Order Date \_\_\_\_\_ Measured By \_\_\_\_\_

Telephone \_\_\_\_\_ Email \_\_\_\_\_ Shipping Method \_\_\_\_\_

Credit Card Info \_\_\_\_\_



**upper leg kit:** Check the box next to the corresponding length and width.

\*cG cannot be more than 30cm larger than cE2.

	Length			Width		Quantity
	Short	Standard	Long	Regular	Wide	
Left	∠ D1-G < 37cm	∠ D1-G 37 - 41cm	∠ D1-G > 41	cG < 80cm	cG < 125cm	
Right				cE2 < 65cm	cE2 < 110cm	

Each upper leg kit includes:

- 1 reduction upper leg component
- 1 Built-In-Tension system guide card
- 1 reduction shelf strap
- 1 paper measuring tape
- 6 hook and loop stays
- 1 direction for use

**knee kit:** (the reduction knee kit is only available in one length and width)

\*cE2 cannot be more than 30cm larger than cC.

	Length	Width	Quantity
Left	30cm	cE2 < 110cm cE < 102cm cC < 80cm	
Right			

Each knee kit includes:

- 1 reduction knee component
- 1 paper measuring tape
- 2 hook and loop knee spines
- 1 direction for use
- 1 reduction shelf strap

**lower leg kit:** Check the box next to the corresponding length and width.

\*cC cannot be more than 30cm larger than cB.

	Length			Width		Quantity
	Short	Standard	Long	Regular	Wide	
Left	∠ mall-D1 < 32cm	∠ mall-D1 32 - 36cm	∠ mall-D1 > 36cm	cC < 65cm	cC < 80cm	
Right				cB < 55cm	cB < 65cm	

Each lower leg kit includes:

- 1 reduction lower leg component
- 1 Built-In-Tension system guide card
- 1 pair reduction kit undersleeve leg
- 1 paper measuring tape
- 1 reduction shelf strap
- 1 direction for use
- 6 hook and loop stays

**foot options:** Indicate quantity in appropriate check boxes.

\*foot options are not included with kits and must be ordered separately.

	customizable interlocking afw	compression anklet		compressive undersocks		pac band™	single band afw
	One Size	Standard	Large	Standard	Large	X-Large	One Size
cA1	< 44cm	19-36 cm	22-41 cm	19-32 cm	24-41 cm	< 44cm	
cB	< 56cm	19-27 cm	25-42 cm	19-27 cm	25-42 cm	< 56cm	
Left							
Right							

**reduction kit extension**

Short (30 cm)	Standard (35 cm)	Long (40 cm)

**additional ancillary items**

x-wide undersleeve leg Quantity (Pairs)	reduction kit shelf strap Quantity (Each)

**foam padding:**

tec pad (2 pack)	flec pad (2 pack)

**reduction kit toe cap\*:** Indicate quantity in appropriate check boxes.

\*toe caps are not included with kits and must be ordered separately.

	Small	Medium	Large	X-Large	XX-Large
cA	20-22 cm	22-24 cm	24-26 cm	26-28 cm	28-30 cm
Left					
Right					

**reduction kit lobe\*:** Indicate quantity in appropriate check boxes.

\*lobe is not included with kits and must be ordered separately.

universal size

Leg Measurement Guide	
G	Groin
E2	15cm above E
E1	5cm above E
E	Center of patella
D1	Knee crease
D	Slightly below knee
C	Widest part of calf
B	Narrowest part of ankle above malleolus
A1	Middle of foot
A	Ball of foot
mall	Center of malleolus
∠	length in a straight line
c	circumference


Customer Name \_\_\_\_\_ Customer No. \_\_\_\_\_ Purchase Order No. \_\_\_\_\_ Fax \_\_\_\_\_

Billing Address \_\_\_\_\_ Patient Name \_\_\_\_\_

Shipping Address \_\_\_\_\_ Order Date \_\_\_\_\_ Measured By \_\_\_\_\_

Telephone \_\_\_\_\_ Email \_\_\_\_\_ Shipping Method \_\_\_\_\_

Credit Card Info \_\_\_\_\_



**arm kit:** Check the box next to the corresponding length and width.

\*If cB is > 30cm hand wrap will not fit, though reduction arm component may still be used.

	Length		Width			Qty.
	Standard	Long	Regular	Wide	X-Wide	
Left	<input type="checkbox"/> C-G < 45cm	<input type="checkbox"/> C-G > 45cm	<input type="checkbox"/> cF < 50cm	<input type="checkbox"/> cF 50-65cm	<input type="checkbox"/> cF 65-80cm	
Right	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> cB < 30cm*	<input type="checkbox"/> cB < 30cm*	<input type="checkbox"/> cB < 30cm*	

Each arm kit includes:

- 1 reduction arm component
- 1 pair undersleeve arm
- 1 customizable hand wrap
- 1 Built-In-Tension guide card
- 1 direction for use

Arm Measurement Guide	
G	Axilla
F	Midpoint of E and G
E	Elbow crease
D	Midpoint of C and E
C	Wrist
B	Base of thumb webbing
ℓ = length in a straight line	
c = circumference	

additional ancillary items:

undersleeve arm standard Standard <50cm max circ. Quantity (Each)	undersleeve arm wide Wide <80cm max circ. Quantity (Each)	reduction kit shelf strap Quantity (Each)
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

reduction kit glove\*: Indicate quantity in appropriate check box.

	Small	Medium	Large	X-Large	XX-Large
cB	18-21 cm	20-23 cm	22-25 cm	24-27 cm	26-29 cm
Left	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Right	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

\*not included in arm kit

Check the box next to the corresponding length and width. Note, width has two measurement criteria. Select the size that meets both criteria.

reduction kit vest:

	Width			Length (Axilla-Waist)		QTY
	Small	Large	Large	Standard	Long	
Chest Circum.	70 - 120cm	121 - 180cm	80 - 180cm	20 - 33cm	34 - 47cm	
Chest Width (Axilla-Axilla Contour)	≤ 39cm	35-39cm	> 39cm			

Each vest kit includes:

- 1 reduction vest component
- 1 flec (flat localized edema control) foam pad
- 1 direction for use

reduction head & neck

Each head & neck kit includes:

- 1 large neck strap
- 1 small facial chin strap
- 1 tec (textured edema control) foam pad
- 1 direction for use

foam padding:

tec pad (2 pack)	flec pad (2 pack)
<input type="checkbox"/>	<input type="checkbox"/>

universal size
<input type="checkbox"/>