



medi

**Silflex**<sup>®</sup>

Soft Silicone Wound  
Healing Solutions

**Adv**ancis  
Medical

medi. I feel better.



# Breaking the Cycle Through Active Healing

## PREPARE

Compress with confidence and prepare wounds with autolytic debridement

# 1

## HEAL

Create optimal wound healing environment and accelerate closure with accurate compression

# 2

## PREVENT

Sustain healthy skin and integrate lifestyle compression to prevent reoccurrence

# 3

## Pain free dressing changes to improve quality of life

- Conforms to the skin's surface for gentle adherence
- Atraumatic dressing designed to minimize pain and trauma
- Large 2mm holes allow for passage of exudate
- Can be cut to size
- Hydrophobic silicone coating prevents dressing from adhering to wound bed



# Silflex®



## Description

Silflex® is an atraumatic wound contact layer made from polyester mesh coated with Silfix® soft silicone.

## Indications

- Acute and chronic wounds that require protections
- Skin Tears
- Skin Abrasions
- Surgical Wounds
- Leg and Pressure ulcers
- Negative Pressure Wound Therapy (as a primary dressing)



Product	SKU	Sizes	Pack Sizes	HCPCS* Code
<b>Silflex®</b> Soft silicone wound contact layer	CR3922	2in x 3in	10	A6206
	CR3923	3in x 4in	10	A6206
	CR3924	5in x 6in	10	A6207
	CR3925	8in x 12in	10	A6208

# Siltape®



## Description

Siltape® is a roll of soft silicone fixation tape (non-sterile) with easy tear perforations, **no scissors are required**.

## Indications

Use in adhesive applications to secure dressings, tubes, IV lines, or for skin tear management. Ideal for pediatrics, newborns and geriatric.

- Scar management, including hypertrophic and keloid scars.
- Protecting tissue from pressure damage or sores from tubing
- Manage skin tears and small incision sites
- Hold dressings in place
- Securing IV lines
- Protect skin from extended use of Personal Protective Equipment (PPE)
- Can be used in any situation where delicate or sensitive skin is an issue



Product	SKU	Sizes	Pack Sizes	HCPCS* Code
<b>Siltape®</b> Soft silicone perforated tape roll	CR3938	.75in x 118in	1	n/a
	CR3939	1.5in x 60in	1	n/a

\*Code usage should be verified by responsible parties in accordance with all published guidelines.



# Silflex<sup>®</sup> Soft Silicone Contact Layer

Pain free dressing changes to improve quality of life

In a recent study Silflex<sup>®</sup> scored better than an alternative silicone contact layer on mean pain scores at application, in situ and removal (see figure 1).<sup>1</sup>



Using Silflex<sup>®</sup> can reduce pain, which improves quality of life, reduces anxiety and may improve the healing process<sup>2</sup>




Figure 1: Mean pain scores

1. Edwards J, Mason S A. Developing a structured process for evaluating burn dressings. Poster presentation  
 2. Holden-Lund C (1987) Effects of relaxation with guided imagery on surgical stress and wound healing. Res Health Nurse 11(4): 235-44

Dressings with Silfix<sup>®</sup> technology can be found in the Eclypse<sup>®</sup>, Advazorb<sup>®</sup> and Silflex<sup>®</sup> ranges.

## Silfix<sup>®</sup> soft silicone dressings



medi USA  
 6481 Franz Warner Parkway  
 Whitsett, NC 27377  
 T 800-633-6334  
 F 888-570-4554  
 mediusa.com

Check out our  
 webinar options!

<https://www.mediusa.com/client-service/education/Webinars-and-Courses/>

Advancis  
 Medical

[www.advancismedical.com](http://www.advancismedical.com)