


# Advazorb

Hydrophilic Foam Dressing Range



Gentle on skin

Distributed by



[www.mediusa.com](http://www.mediusa.com)

**Adv**  **ancis**  
Medical

[www.advancismedical.com](http://www.advancismedical.com)

# Advazorb®

## Hydrophilic foam dressing range

Manufactured in the UK, the Advancis Medical® range of Advazorb® Hydrophilic foam dressings, offer the combination of a high quality performance product at a competitive price.

Advazorb® is a comprehensive range of patient friendly, absorbent foam dressings presented in non-adhesive and atraumatic silicone adhesive format. The range is designed to overcome the complex challenges of managing exudate whilst protecting 'at risk' fragile skin, a combination of excellent fluid handling capabilities of Advazorb® with a layer of atraumatic soft silicone adhesive.

**Advazorb®, Advazorb Silfix® and Advazorb Border® are available in both a regular thickness for moderate to high exudate and 'Lite' versions for low to moderately exuding wounds.**



**Advazorb®**  
Hydrophilic foam dressing  
with film backing

**Soft and conformable, non-adherent, hydrophilic, polyurethane foam dressings with a breathable film backing.**

A non-adherent foam dressing with excellent fluid handling properties. Ideal for use under compression bandaging with its fluid retention, low profile and low friction film backing.

- High quality polyurethane foam, rapidly absorbs and retains what it absorbs. Available in regular and 'Lite' versions to suit a wider range of wounds.
- The Advazorb® range features a waterproof film backing that prevents strike-through. The backing has been designed to reduce friction on clothes and bandages which extends dressing wear time by reducing 'rucking' and preventing dressings falling off.
- With an MVTR range of 2600-2700 g/m<sup>2</sup>/24hours the Advazorb® foam range provides the ideal conditions for wound healing.
- Increasing patient compliance and comfort.
- Can be used on a wide variety of exuding wounds including under compression and can be cut to shape and size.



**Advazorb Heel®**  
Hydrophilic foam dressing  
with film backing

**A conformable, non-adherent, hydrophilic, polyurethane foam dressing with a breathable film backing designed to fit around the heel.**

A non-adherent foam dressing with excellent fluid handling properties. Ideal for use under compression bandaging with its fluid retention, low profile and low friction film backing.

- The same high quality PU foam as the other dressings in the Advazorb® range, it has excellent fluid handling properties and rapidly absorbs and retains exudate.
- It features the low friction, bacteria-proof, thin film backing that the entire Advazorb® range has meaning extended wear time, reduced 'rucking' and is waterproof.
- Soft and conformable, increasing patient comfort and is ideal to use under compression due to its low profile.
- It is anatomically shaped to fit around the heel so that the dressing is in close contact with the wound to provide an environment conducive to healing.

Product	Size	Dressings per box	HCPCS code
Advazorb®	3in x 3in	10	A6209
	4in x 4in	10	A6209
	5in x 5in	10	A6210
	6in x 6in	10	A6210
	4in x 8in	10	A6210
	8in x 8in	10	A6211
Advazorb Heel®	6.7in x 8.3in	5	A6210

**medi** USA  
6481 Franz Warner Pkwy  
27377 Whitsett  
United States of America  
Tel: 800-633-6334 | E-mail: info@mediusa.com

**Advancis**  
Medical  
www.advancismedical.com