


NEW

# Advazorb

## Areola

### The protective and secure breast dressing developed with clinicians

- No more 'patchwork' dressings
- Low friction backing to prevent uncomfortable snagging on clothes
- Reduced dressing time and associated patient discomfort
- Central hole to prevent pressure on the nipple including incisions to accommodate larger nipples
- Designed for patient comfort
- Silfix soft silicone adhesive for atraumatic removal




# An Innovative Breast Dressing

**Advazorb Areola** is the result of a collaboration with Nottingham University Hospitals and Nottingham Breast Institute to fulfil a very specific need – a protective and secure breast dressing.

**Advazorb Areola** combines Silfix® soft silicone with conformable, absorbent foam to deliver excellent patient comfort, effective exudate management, secure fixation and atraumatic removal.


## Three layers of performance

- 1 Low friction backing**  
Increases 'wearability' and prevents uncomfortable snagging on clothes
- 2 Silfix silicone adhesive**  
For easy placement and replacement  
Atraumatic removal
- 3 Comfortable, absorbent foam**  
With nipple hole to prevent nipple pressure.  
Incisions ensure the hole suits larger nipple sizes.



## The challenge of dressing the breast

Clinicians are often faced with trying to create their own breast dressing from a patchwork of square dressings that do not conform to the shape of the breast.



### Fit for purpose

**Advazorb Areola** is ergonomically designed to fit the contours of the breast making dressing application, easy, quick and comfortable for the patient.

## The ideal choice for:

- Mastopexy
- Breast reduction
- Nipple reconstruction
- Breast augment through areola incision
- Gynaecomastia
- Other Nipple / Areola surgery



## Dimensions:

Total width **5 in.**

Foam diameter **1.96 in.**

Nipple hole **0.5 in. plus incisions**

## Product Code:

**CR4490**

For further information or to request a demonstration on the innovative Advazorb Areola, contact medi USA at **800-633-6334** or email **WoundCare@mediusa.com**



medi USA  
6481 Franz Warner Parkway  
Whitsett, NC 27377  
T 800-633-6334  
F 888-570-4554  
mediusa.com