

BC-55 Upgrade, Restore, Detection Level Adj

1. BC55 software upgrade

1.1 download & copy software to USB

Download & Copy the software (a file with extension .ZPK) to USB memory stick.

For example, copy a file "BC55_06010301046_06_20211213.E0644173534F5AEB.ZPK" to the USB **root** directory,

Note: don't change the file name and make sure there is only one .ZPK software under this directory.

1.2 upgrade steps

1) turn off the machine, plug USB, then turn on the machine



2) the machine will recognize the software and upgrade itself automatically as below



Please switch off the machine, then take away USB, and reboot it.